

# Mouse Esm1/ Endothelial cell-specific molecule 1 ELISA Kit

Catalog No: #EK21129

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK21129-1 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Esm1/ Endothelial cell-specific molecule 1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Uniprot	Q9QYY7

## Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:8.02pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only