

## Recombinant Human Zinc finger protein ZPR1 (ZPR1)

Catalog No: #AP89685



Package Size: #AP89685-1 20ug #AP89685-2 100ug #AP89685-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Zinc finger protein ZPR1 (ZPR1)
Brief Description	Recombinant Proteins
Uniprot	O75312
Target Species	Homo sapiens (Human)
SDS-PAGE MW	58.4 kDa
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Target Sequence	<p>Protein Length: Full length</p> <p>Expression Region: 1-459aa</p> <p>MAASGAVEPGPPGAAVAPSPAPAPPPAPDHLFRPISAEDDEEQPTEIESLCMNCYCNMTRLLLLTKIPFFREII  VSSFSCHEHCGWNNTEIQSAGRIQDQGVRYTLSVRALEDMNREVVKTDSAATRIPELDFEIPAFSQKGALTTVE  GLITRAISGLEQDQPARRANKDATAERIDEFIVKLKELKQVASPFTLIIDDPGNSFVENPHAPQKDDALVITHYN  RTRQQEEMGLQEEAPAEKPEEEDLRNEVLQFSTNCPECNAPAQTNMMLVQIPHFKVEIIMATNCENGHRT  NEVKSGGAVEPLGTRITLHITDASDMTRDLLKSETCSVEIPELEFELGMAVLGGKFTTLEGLLDIRELVTKNPF  TLGDSSNPGQTERLQEFQKMDQIEGNMKAHFIMDDPAGNSYLQNVYAPEDDPKVERYKRTFDQNEELG  LNDMKTEGYEAGLAPQR</p>
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

## Application Details

Greater than 90% as determined by SDS-PAGE.

Note: This product is for in vitro research use only