

Arrestin- $\beta$ -1 Polyclonal Antibody

Catalog No: #48857

Package Size: #48857-1 50ul #48857-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Arrestin- $\beta$ -1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Arrestin- $\beta$ -1 Polyclonal Antibody detects endogenous levels of Arrestin- $\beta$ -1 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human ARRB1.
Target Name	ARRB1
Other Names	ARRB1; ARR1; Beta-arrestin-1; Arrestin beta-1
Uniprot	P49407
GeneID	408
Calculated MW	50kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1:500 - 1:2000. IHC 1:100 - 1:300. IF 1:200 - 1:1000. ELISA: 1:10000.

## Images



Western blot analysis of lysates from Jurkat and HUVEC cells, using ARRB1 Antibody. The lane on the right is blocked with the synthesized peptide.

---

Note: This product is for in vitro research use only