

## TRAP Mouse mAb

Catalog No: #60450

Package Size: #60450-1 50ul #60450-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	TRAP Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Applications	IHC;IF
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody can specifically recognize human TRAP protein.
Immunogen Description	Synthesized peptide derived from human TRAP
Target Name	ACP5
Other Names	Tartrate-resistant acid phosphatase type 5 (TR-AP;EC 3.1.3.2;Tartrate-resistant acid ATPase;TrATPase;Type 5 acid phosphatase)
Uniprot	P13686
GeneID	54
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

IHC-p 1:200-400, IF 1:200-400

## Images



Human lung tissue was stained with Anti-TRAP Antibody



Human giant cell tumor tissue was stained with Anti-TRAP Antibody

---

Note: This product is for in vitro research use only