

MT1A rabbit pAb

Catalog No: #48877

Package Size: #48877-1 50ul #48877-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MT1A rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC;IF
Species Reactivity	Human;Rat;Mouse
Specificity	This antibody detects endogenous levels of MT1A at Human
Immunogen Description	Synthesized peptide derived from human MT1A
Target Name	MT1A MT1S
Uniprot	P04731
GeneID	4489
Calculated MW	7kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

IHC-p 1 μ O 50-200. IF 1:50-200

Images



Immunohistochemical analysis of paraffin-embedded human liver cancer.

Note: This product is for in vitro research use only