

Mouse Pf4/ Platelet factor 4 ELISA Kit

Catalog No: #EK21120



Package Size: #EK21120-1 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Pf4/ Platelet factor 4 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Other Names	C-X-C motif chemokine 4, PF-4, Cxcl4, Scyb4
Uniprot	Q9Z126

Application Details

Detect Range:78-5000pg/mL

Sensitivity:39.8pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only