

Porcine Caspase 3 (CASP3) ELISA Kit

Catalog No: #EK20613



Package Size: #EK20613-1 48T #EK20613-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|-------------------------------------|
| Product Name | Porcine Caspase 3 (CASP3) ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Porcine |

Application Details

Detect Range:0.312-20 ng/mL

Sensitivity:0.105ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only