

SELENOP Rabbit pAb

Catalog No: #55654

Package Size: #55654-1 50ul #55654-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SELENOP Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB IF/ICC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human SELENOP.
Other Names	SeP; SELP; SEPP; SEPP1; SELENOP
Uniprot	P49908
GeneID	6414
Calculated MW	43KDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:1000 IF/ICC 1:100 - 1:500

Images



Western blot analysis of Rat liver, using SELENOP Rabbit pAb at 1:1000 dilution.



Immunofluorescence analysis of NIH/3T3 cells using SELENOP Rabbit pAb at dilution of 1:500 (40x lens). Blue: DAPI for nuclear staining.

Note: This product is for in vitro research use only