

FUSIP1 rabbit pAb

Catalog No: #48878



Package Size: #48878-1 50ul #48878-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FUSIP1 rabbit pAb
Brief Description	rabbit polyclonal antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse
Specificity	This antibody detects endogenous levels of FUSIP1 at Human, Mouse
Immunogen Description	Synthesized peptide derived from human FUSIP1
Target Name	SRSF10
Other Names	SRSF10 FUSIP1 FUSIP2 SFRS13A TASR
Uniprot	O75494
GeneID	100996657
Calculated MW	29kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only