

FRPD2 rabbit pAb

Catalog No: #48886

Package Size: #48886-1 50ul #48886-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FRPD2 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Rat;Mouse
Specificity	This antibody detects endogenous levels of FRPD2 at Human
Immunogen Description	Synthesized peptide derived from human FRPD2
Target Name	FRMPD2 PDZD5C PDZK4 PDZK5C
Uniprot	Q68DX3
GeneID	143162
Calculated MW	144kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1αO 500-2000

Images



Western blot analysis of lysates from THP-1 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only