

HNF1 α (Phospho Ser247) rabbit pAb

Catalog No: #14303



Package Size: #14303-1 50ul #14303-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HNF1 α (Phospho Ser247) rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human;Rat;Mouse
Specificity	This antibody detects endogenous levels of Human HNF1 α (phospho-Ser247)
Immunogen Description	Synthesized phospho peptide around human HNF1 α (Ser247)
Target Name	HNF1A TCF1
Other Names	Hepatocyte nuclear factor 1-alpha (HNF-1-alpha) (HNF-1A) (Liver-specific transcription factor LF-B1) (LFB1) (Transcription factor 1) (TCF-1)
Uniprot	P20823
GeneID	6927
Calculated MW	69kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only