

DAPK1 (Phospho Ser308) Polyclonal Antibody

Catalog No: #14304

Package Size: #14304-1 50ul #14304-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DAPK1 (Phospho Ser308) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC;IF;WB
Species Reactivity	Human;Mouse
Specificity	This antibody detects endogenous phospho levels of DAPK1 (Phospho-Ser308) at Human:S308, Mouse:S308
Immunogen Description	Synthesized peptide derived from human DAPK1 (Phospho-Ser308)
Target Name	DAPK1 DAPK
Other Names	Death-associated protein kinase 1 (DAP kinase 1;EC 2.7.11.1)
Uniprot	P53355
GeneID	1612
Calculated MW	160kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

IHC-p 1:50-200, WB 1:500-2000. IF 1:50-200

Images



Immunohistochemical analysis of paraffin-embedded human tonsil.

Note: This product is for in vitro research use only