

# Polycystin 1 rabbit pAb

Catalog No: #48915

Package Size: #48915-1 50ul #48915-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Polycystin 1 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse
Specificity	This antibody detects endogenous levels of Polycystin 1 at Human, Mouse
Immunogen Description	Synthesized peptide derived from human Polycystin 1
Target Name	PKD1
Uniprot	P98161
GeneID	5310
Calculated MW	473kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1:500-2000

Note: This product is for in vitro research use only