

UBE3C Polyclonal Antibody

Catalog No: #48924

Package Size: #48924-1 50ul #48924-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	UBE3C Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;ELISA
Species Reactivity	Human;Rat;Mouse
Specificity	UBE3C Polyclonal Antibody detects endogenous levels of UBE3C protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human UBE3C.
Target Name	UBE3C
Other Names	UBE3C; KIAA0010; KIAA10; Ubiquitin-protein ligase E3C; HectH2
Uniprot	Q15386
GeneID	148581
Calculated MW	125kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500 - 1:2000. ELISA: 1:10000

Images



Western Blot analysis of HuvEc cells using UBE3C Polyclonal Antibody.

Note: This product is for in vitro research use only