

CD204 Polyclonal Antibody

Catalog No: #48940



Package Size: #48940-1 50ul #48940-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD204 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous CD204
Immunogen Description	The antiserum was produced against synthesized peptide derived from the Internal region of human MSR1.
Target Name	MSR1 SCARA1
Other Names	Macrophage scavenger receptor types I and II (Macrophage acetylated LDL receptor I and II;Scavenger receptor class A member 1;CD antigen CD204)
Uniprot	P21757
GeneID	4481
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

IHC-p 1:50-200, ELISA 1:10000-20000. IF 1:50-200

Note: This product is for in vitro research use only