

Mouse Glia-derived nexin(SERPINE2) ELISA kit

Catalog No: #EK20628



Package Size: #EK20628-1 48T #EK20628-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Glia-derived nexin(SERPINE2) ELISA kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.312-20ng/ml

Sensitivity:0.112ng/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only