

PIK3CA Phospho Tyr317 rabbit pAb

Catalog No: #14305



Package Size: #14305-1 50ul #14305-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PIK3CA Phospho Tyr317 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of PIK3CA Phospho-Tyr317 at Human, Mouse,Rat
Immunogen Description	Synthesized peptide derived from human PIK3CA Phospho-Tyr317
Target Name	PIK3CA
Uniprot	P42336
GeneID	5290
Calculated MW	117kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only