

## LATS1/2 (Phospho Thr1079/1041) Antibody

Catalog No: #14306



Package Size: #14306-1 50ul #14306-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	LATS1/2 (Phospho Thr1079/1041) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	This detects endogenous levels of LATS1/2 (Phospho-Thr1079/1041)
Immunogen Description	Synthesized phospho derived from human LATS1/2 (Phospho-Thr1079/1041)
Target Name	LATS1 WARTS
Other Names	Serine/threonine-protein kinase LATS1 (EC 2.7.11.1) (Large tumor suppressor homolog 1) (WARTS protein kinase) (h-warts)
Uniprot	O95835
GeneID	9113
Calculated MW	140kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1:500-2000, ELISA 1:10000-20000

Note: This product is for in vitro research use only