

ABCA1 Mouse mAb

Catalog No: #52906



Package Size: #52906-1 50ul #52906-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ABCA1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Purification	Affinity Purified
Applications	ICC/IF,FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human ABCA1 expressed in E. Coli.
Target Name	ABCA1
Other Names	TGD; ABC1; CERP; ABC-1; HDLDT1; HPALP1; HDLCQTL13
Uniprot	O95477
GeneID	19
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IF: 1/50-1/200 FC: 1/50-1/100

Note: This product is for in vitro research use only