

# Mouse Membrane-bound transcription factor site-1 protease (MBTPS1) ELISA Kit



Catalog No: #EK20641

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20641-1 48T #EK20641-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Membrane-bound transcription factor site-1 protease (MBTPS1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

## Application Details

Detect Range:0.312-20ng/mL
Sensitivity:0.108ng/ml
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 &mu;L
Assay Time:1-4.5h
Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only