

SP130 Polyclonal Antibody

Catalog No: #48945



Package Size: #48945-1 50ul #48945-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SP130 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;ELISA
Species Reactivity	Human;Mouse
Specificity	SP130 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Description	Synthesized peptide derived from part region of human protein
Target Name	SAP130
Uniprot	Q9H0E3
GeneID	79595
Calculated MW	115kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000 ELISA 1:5000-20000

Note: This product is for in vitro research use only