

ZNF580 Rabbit pAb

Catalog No: #48947



Package Size: #48947-1 50ul #48947-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ZNF580 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	affinity purified by Protein A
Applications	WB, IHC-P, IHC-F, ICC, IF
Species Reactivity	Human;Mouse;Rat
Immunogen Description	KLH conjugated synthetic peptide derived from human ZNF580: 51-150/172
Conjugates	Unconjugated
Other Names	LDL induced EC protein; LDL-induced EC protein; Zinc finger protein 580; ZN580_HUMAN; Znf580.
Uniprot	Q9UK33
GeneID	51157
Concentration	1mg/ml
Formulation	Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000

IHC-P 1:100-500

IHC-F 1:100-500

ICC 1:100-500

IF 1:100-500

Note: This product is for in vitro research use only