

GPHB5 Rabbit pAb

Catalog No: #48970



Package Size: #48970-1 50ul #48970-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GPHB5 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	affinity purified by Protein A
Applications	IHC-P, IHC-F, ICC, IF, ELISA
Species Reactivity	Human;Mouse;Rat
Immunogen Description	KLH conjugated synthetic peptide derived from human GPHB5: 25-100/130
Other Names	CGH; Glycoprotein hormone beta 5; GPHB5_HUMAN; Glycoprotein hormone beta 5 precursor; Glycoprotein hormone beta subunit; GPB5; Thyrostimulin subunit beta; ZLUT1.
Uniprot	Q86YW7
GeneID	122876
Calculated MW	12kD
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P 1:100-500

IHC-F 1:100-500

ICC 1:100-500

IF 1:100-500

ELISA 1:5000-10000

Note: This product is for in vitro research use only