

ALDH3A1 Rabbit pAb

Catalog No: #49026



Package Size: #49026-1 50ul #49026-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | ALDH3A1 Rabbit pAb |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | affinity purified by Protein A |
| Applications | WB, IHC-P, IHC-F, ICC, IF, ELISA |
| Species Reactivity | Human;Mouse;Rat |
| Immunogen Description | KLH conjugated synthetic peptide derived from human ALDH3A1: 151-250/453 |
| Other Names | AHD 4; AHD C; Ahd4; AHDC; AL3A1_HUMAN; Aldehyde dehydrogenase 3; Aldehyde dehydrogenase 3 family member A1; Aldehyde dehydrogenase 3 family, member A1; Aldehyde dehydrogenase 3A1; Aldehyde dehydrogenase; Aldehyde dehydrogenase class 3; Aldehyde dehydrogenase dimeric NADP preferring; Aldehyde dehydrogenase family 3 member A1; Aldehyde dehydrogenase family 3 subfamily A1; Aldehyde dehydrogenase isozyme 3; Aldehyde dehydrogenase type III; Aldehyde dehydrogenase, dimeric NADP preferring; ALDH 3; ALDH 3A1; ALDH III; ALDH3 A1; ALDH3; ALDH3A 1; Aldh3a1; ALDHIII; dimeric NADP-preferring; MGC10406; Stomach aldehyde dehydrogenase; Tumor associated aldehyde dehydrogenase tumor ALDH. |
| Uniprot | P30838 |
| GeneID | 218 |
| Calculated MW | 50kD |
| Concentration | 1mg/ml |
| Formulation | 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB 1:500-2000

IHC-P 1:100-500

IHC-F 1:100-500

ICC 1:100-500

IF 1:100-500

ELISA 1:5000-10000

Note: This product is for in vitro research use only