

Fish Thyroid stimulating hormone (TSH) ELISA Kit

Catalog No: #EK20652



Package Size: #EK20652-1 48T #EK20652-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fish Thyroid stimulating hormone (TSH) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Fish

Application Details

Detect Range: 0.3-12 μ g/mL

Sensitivity: 0.1 μ g/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only