

Recombinant Mouse Galectin-1 (Lgals1)

Catalog No: #AP89677



Package Size: #AP89677-1 20ug #AP89677-2 100ug #AP89677-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Galectin-1 (Lgals1)
Brief Description	Recombinant Proteins
Host Species	E.coli
Uniprot	P16045
Target Species	Mus musculus (Mouse)
Calculated MW	22.2 kDa
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Target Sequence	Protein Length:Full Length of Mature Protein Expression Region: 2-135aa ACGLVASNLNLKPGECKVRGEVASDAKSFVLNLGKDSNNLCLHFNPRFNAHGDANTIVCNTKEDGTWGTEH REPAFPFQPGSITEVCITFDQADLTIKLPDGHEFKFPNRLNMEAINYMAADGDFKIKCVAFE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Application Details

Greater than 90% as determined by SDS-PAGE.

Note: This product is for in vitro research use only