

DTX4 rabbit pAb

Catalog No: #49028

Package Size: #49028-1 50ul #49028-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DTX4 rabbit pAb
Brief Description	rabbit polyclonal antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse
Specificity	This antibody detects endogenous levels of DTX4 at Human, Mouse
Immunogen Description	Synthesized peptide derived from human DTX4
Target Name	DTX4 KIAA0937 RNF155
Uniprot	Q9Y2E6
GeneID	23220
Calculated MW	68kD
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only