

# Rat SH3 and multiple ankyrin repeat domains protein 3 (SHANK3) ELISA Kit

Catalog No: #EK20661

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20661-1 48T #EK20661-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat SH3 and multiple ankyrin repeat domains protein 3 (SHANK3) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat

## Application Details

Detect Range:62.5-4000pg/mL

Sensitivity:23.4pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50  $\mu$ L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only