

Recombinant Human Developmental pluripotency-associated protein 2 (DPPA2)

Catalog No: #AP89679

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP89679-1 20ug #AP89679-2 100ug #AP89679-3 1mg

Description

| | |
|-------------------|---|
| Product Name | Recombinant Human Developmental pluripotency-associated protein 2 (DPPA2) |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Accession No. | Q7Z7J5 |
| Target Species | Homo sapiens (Human) |
| Target Sequence | <p>Expression Region:1-298aa</p> <p>Sequence Info:Full Length</p> <p>MSDANLDSSK KNFLEGEVDD EESVILTLVP VKDDANMEQM EPSVSSTSDV KLEKPKKYNP GHLLQTNEQF TAPQKARCKI PALPLPTILP PINKVCRDTL RDWCQQLGLS TNGKKIEVYL RLHRHAYPEQ RQDMPEMSQE TRLQRCRKR KAVTKRARLQ RSYEMNERAE ETNTVEVITS APGAMLASWA RIAARAVQPK ALNSCSIPVS VEAFMQASG VRWCVVHGRL LSADTKGWVR LQFHAGQAVV PTTHRRMISL FLLPACIFPS PGIEDNMLCP DCAKRNKMM KRLMTVEK</p> |
| Storage | Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles. |

Application Details

Greater than 85% as determined by SDS-PAGE.

Note: This product is for in vitro research use only