

MPL Antibody FITC Conjugated

Catalog No: #C00717F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MPL Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC,IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 410-445 635 derived from human MPL
Conjugates	FITC
Target Name	MPL
Other Names	MPLV; TPOR; C-MPL; CD110; THCYT2; Thrombopoietin receptor; TPO-R; Myeloproliferative leukemia protein; Proto-oncogene c-Mpl; MPL
Accession No.	Swiss-Prot#P40238NCBI Gene ID4352
Uniprot	P40238
GeneID	4352;
Excitation Emission	494nm 518nm
Cell Localization	Helical
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

Receptor for thrombopoietin. May represent a regulatory molecule specific for TPO-R-dependent immune responses.

Note: This product is for in vitro research use only