

Rat Atf6/ Cyclic AMP-dependent transcription factor ATF-6 alpha ELISA Kit

Catalog No: #EK21130

Orders: order@signalwayantibody.com

Package Size: #EK21130-1 96T

Support: tech@signalwayantibody.com

Description

Product Name	Rat Atf6/ Cyclic AMP-dependent transcription factor ATF-6 alpha ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat
Other Names	Activating transcription factor 6 alpha, ATF6-alpha, cAMP-dependent transcription factor ATF-6 alpha
Uniprot	G3V909

Application Details

Detect Range:0.156-10ng/mL
Sensitivity:0.074ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only