

Rat Eif2ak3/ Eukaryotic translation initiation factor 2-alpha kinase 3 ELISA Kit



Catalog No: #EK21135

Orders: order@signalwayantibody.com

Package Size: #EK21135-1 96T

Support: tech@signalwayantibody.com

Description

Product Name	Rat Eif2ak3/ Eukaryotic translation initiation factor 2-alpha kinase 3 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat
Other Names	PRKR-like endoplasmic reticulum kinase, Pancreatic eIF2-alpha kinase, Pek, Perk
Uniprot	Q9Z1Z1

Application Details

Detect Range:0.156-10ng/mL
Sensitivity:0.081ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only