

Rat Endoplasmic Reticulum To Nucleus Signaling 1 (ERN1) ELISA Kit

Catalog No: #EK20657

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20657-1 48T #EK20657-2 96T

Description

Product Name	Rat Endoplasmic Reticulum To Nucleus Signaling 1 (ERN1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:62.5-4000 pg/mL

Sensitivity:25.7pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only