

Human Ig alpha-1 chain C region (IGHA1) ELISA Kit

Catalog No: #EK20653



Package Size: #EK20653-1 48T #EK20653-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Ig alpha-1 chain C region (IGHA1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range: 1.56-100 µg/mL

Sensitivity: 0.52 µg/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 µL

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only