

Human Bestrophin-1 (BEST1) ELISA Kit

Catalog No: #EK20675



Package Size: #EK20675-1 48T #EK20675-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Bestrophin-1 (BEST1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:125-8000pg/mL

Sensitivity:45pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture

supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only