

Mouse Protein phosphatase 1 regulatory subunit 15A (PPP1R15A) ELISA Kit



Catalog No: #EK20669

Orders: order@signalwayantibody.com

Package Size: #EK20669-1 48T #EK20669-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Protein phosphatase 1 regulatory subunit 15A (PPP1R15A) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:31.2-2000 pg/mL

Sensitivity:12.5pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only