

# Rat Protein phosphatase 1 regulatory subunit 15A (PPP1R15A) ELISA Kit



Catalog No: #EK20670

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20670-1 48T #EK20670-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Protein phosphatase 1 regulatory subunit 15A (PPP1R15A) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat

## Application Details

Detect Range:15.6-1000 pg/mL
Sensitivity:4.75pg/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 &mu;L
Assay Time:1-4.5h
Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only