

GBG11 rabbit pAb

Catalog No: #49064

Package Size: #49064-1 50ul #49064-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GBG11 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC;IF
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of GBG11 at Human/Mouse/Rat
Immunogen Description	Synthesized peptide derived from human GBG11
Target Name	GNG11 GNGT11
Uniprot	P61952
GeneID	2791
Calculated MW	8kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

IHC-p 1 μ O 50-200. IF 1:50-200

Images



Immunohistochemical analysis of paraffin-embedded human tonsil.

Note: This product is for in vitro research use only