

## SCNM1 Polyclonal Antibody

Catalog No: #49083

Package Size: #49083-1 50ul #49083-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	SCNM1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	SCNM1 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Description	Synthesized peptide derived from part region of human protein
Target Name	SCNM1
Uniprot	Q9BWG6
GeneID	100534012
Calculated MW	25kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1:500-2000 ELISA 1:5000-20000

## Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only