

Sheep Interleukin 5 (IL5) ELISA Kit

Catalog No: #EK20685



Package Size: #EK20685-1 48T #EK20685-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sheep Interleukin 5 (IL5) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Sheep

Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:7.8pg/mL

Sample Type:serum, plasma, tissue homogenates or other biological fluids.

Sample Volume: 100μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only