

# Human 45 kDa calcium-binding protein (SDF4) ELISA Kit



Catalog No: #EK20690

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20690-1 48T #EK20690-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human 45 kDa calcium-binding protein (SDF4) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

## Application Details

Detect Range:0.156-10 ng/mL

Sensitivity:0.067ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only