

# Canine C-Reactive Protein (CRP) ELISA Kit

Catalog No: #EK20697



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Canine C-Reactive Protein (CRP) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Canine

## Application Details

Detect Range: 0.312-20  $\mu$ g/mL

Sensitivity: 0.105  $\mu$ g/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50  $\mu$ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

## Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only