

Sarcomeric Alpha Actinin Antibody FITC Conjugated

Catalog No: #C00719F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sarcomeric Alpha Actinin Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt, ICC, IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 480-530 894 derived from human Sarcomeric Alpha Actinin
Conjugates	FITC
Target Name	Sarcomeric Alpha Actinin
Other Names	CMD1AA; Alpha-actinin-2; Alpha-actinin skeletal muscle isoform 2; F-actin cross-linking protein; ACTN2
Accession No.	Swiss-Prot#P35609NCBI Gene ID88
Uniprot	P35609
GeneID	88;
Excitation Emission	494nm 518nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test ICC=1:50-200 IF=1:50-200

Background

F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling protein.

Note: This product is for in vitro research use only