

Rat Rheumatoid Factor (RF) ELISA Kit

Catalog No: #EK20699



Package Size: #EK20699-1 48T #EK20699-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--------------------------------------|
| Product Name | Rat Rheumatoid Factor (RF) ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Rat |

Application Details

Detect Range:250-4000 pg/mL

Sensitivity:125pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only