

# Rat Anti-cyclic Citrullinated Peptide Antibody (CCP-Ab) ELISA Kit



Catalog No: #EK20703

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20703-1 48T #EK20703-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Anti-cyclic Citrullinated Peptide Antibody (CCP-Ab) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat

## Application Details

Detect Range:3.12-200 IU/mL

Sensitivity:1.27 IU/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only