

Human Choline-phosphate cytidyltransferase A (PCYT1A) ELISA Kit



Catalog No: #EK20708

Orders: order@signalwayantibody.com

Package Size: #EK20708-1 48T #EK20708-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Human Choline-phosphate cytidyltransferase A (PCYT1A) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20ng/mL
Sensitivity:0.126ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only