

# Human CMP-N-acetylneuraminase-poly-alpha-2,8-sialyltransferase (ST8SIA4) ELISA Kit

Catalog No: #EK20716

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #EK20716-1 48T #EK20716-2 96T

## Description

Product Name	Human CMP-N-acetylneuraminase-poly-alpha-2,8-sialyltransferase (ST8SIA4) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

## Application Details

Detect Range:0.156-10ng/mL
Sensitivity:0.062ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 &mu;L
Assay Time:1-4.5h
Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only