

## VEGFA Mouse mAb

Catalog No: #49212

Package Size: #49212-1 50ul #49212-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	VEGFA Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Applications	WB;ELISA
Species Reactivity	Human
Specificity	VEGF-A Antibody detects endogenous levels of human VEGF-A protein.. The antibody was also Predict react with Mouse;Rat
Immunogen Description	The antiserum was produced against synthesized peptide derived from human VEGF-A.
Target Name	VEGFA
Other Names	VEGFA;VEGF;Vascular endothelial growth factor A;VEGF-A;Vascular permeability factor;VPF
Uniprot	P15692
GeneID	7422
Calculated MW	21kD(monomer),42kD(dimer)
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1:500-2000,ELISA 1:5000-20000

## Images



Whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with anti-VEGF-A antibody.

---

Note: This product is for in vitro research use only