

Porcine Secretory immunoglobulin A (SIgA) ELISA Kit

Catalog No: #EK20715



Package Size: #EK20715-1 48T #EK20715-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Porcine Secretory immunoglobulin A (SIgA) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Porcine

Application Details

Detect Range: 1.56-100 μ g/mL

Sensitivity: 0.51 μ g/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only