

Human Cross Linked C-Telopeptide Of Type I Collagen (CTXI) ELISA Kit



Catalog No: #EK20546

Orders: order@signalwayantibody.com

Package Size: #EK20546-1 48T #EK20546-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Human Cross Linked C-Telopeptide Of Type I Collagen (CTXI) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.156-10ng/mL
Sensitivity:0.059ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm
Detection Method:sandwich

Note: This product is for in vitro research use only